




# BIOREACTOR

Bioreaction Engineering



Pharmabios®

E: [info@pharmabios.com](mailto:info@pharmabios.com)

W: [www.pharmabios.com](http://www.pharmabios.com)

## COMPANY PRESENTATION / BIOREACTORS

---

**Pharmabios**<sup>®</sup> is a company focused on the design and manufacturing of high technology equipment for the pharmaceutical and biotech branch. It is also a supplier for engineering services on these fields.

The **Pharmabios**<sup>®</sup> bioreactors are specially designed for the scale-up and manufacturing of biological pharmaceutical products. They are equipped with all the necessary components to carry out the bioreaction.



## WORKING PRINCIPLE:

---

The **Pharmabios**<sup>®</sup> bioreactors are designed to work in “batch” mode, with the following operational steps:

**Load of the culture media:** The required media volume is added to the bioreactor.


**Inoculum addition:** The inoculum is added to the bioreactor.

**Bioreaction:** The bioreaction is performed, where the temperature, optical density, pH and dissolved oxygen are controlled.

**Bioreactor unloading:** The product is unloaded to process it in the downstream steps.

**CIP:** Automated cleaning of the bioreactor and its piping and components.

**SIP:** Bioreactor sterilisation.




## PHARMABIOS® EQUIPMENT:

**General features:** Finish of the product-contact parts  $Ra \leq 0.5 \mu m.$ , mechanically polished. Valves in product-contact are diaphragm valves with diaphragm in EPDM. The design minimises the use of connections and welds. Components in product-contact suitable for the pharmaceutical industry. Design carried out in order to minimize the number of connections and welds, dead legs comply with the standard  $L < 3D$ , product contact components designed for the pharmaceutical industry.

**Supporting frame:** All the components including the bioreactor are supported on a stainless steel 304L frame.

**Control system:** The equipment is controlled by means of an electrical cabinet which contains all the necessary components for the equipment performance. The control system is composed of a touch screen HMI (SCADA) type and a Siemens PLC.

**Included documentation:** User and maintenance manuals, quality certificates of the components, welding certificates, manuals/data sheets of the components, CE marking DS/DQ/FAT protocols to be executed in our facilities.




In PHARMABIOS® we adapt ourselves to the customer needs, offering custom-made solutions.

## THE ADVANTAGES OF THE PHARMABIOS® DESIGN:

- » A complete quality documentation, including FDS protocols, DQ & FAT (IQ/OQ) as option.
- » A very versatile control system that allows the configuration of any options to carry out the culture in optimal conditions.
- » A high quality equipment, with GMP design and compliant with all the current international normatives (21 CFR Part 11, GAMP 5, ASME BPE and ISPE recommendations).
- » Customised technical consultancy.
- » Possibility of customisation, tailor-made.
- » High range components, with recognised brands.
- » Sanitary connections, ASME BPE.

## OPTIONS:

- » Heating system using resistances.
- » Setup for anaerobic growth (N<sub>2</sub> addition, CO<sub>2</sub>, ...).
- » Control system upgrade.
- » Load cells.
- » Interface with other equipment or components.
- » Functioning of the bioreactor in continuous mode.
- » Extra copies for documentation.
- » 6-channel graphic register.
- » Supplementary instruments.



## TECHNICAL SPECIFICATIONS:

MODEL	BIO-50	BIO-100	BIO-500	BIO-1000	BIO-2000
Bioreactor volume (l)	50	100	500	1000	2000
Dimensions (length x height x width)	1.500 x 1.500 x 2.000	2.000 x 2.000 x 2.500	2.500 x 2.500 x 2.500	3.000 x 2.500 x 3.000	3.500 x 3.000 x 3.500
PLC/Touch screen	Siemens S7-1200, Siemens, con Tia-portal				
<b>Utilities required</b>					
Electrical consumption (kW)	2	3	4	6	11
Clean steam (flow rate / pressure)	15 kg/h / 3-8bg	30 kg/h / 3-8bg	50 kg/h / 3-8bg	80 kg/h / 3-8bg	100 kg/h / 3-8bg
Purified water (flow rate/pressure)	0.5 m <sup>3</sup> /h / 2bg	1 m <sup>3</sup> /h / 2bg	1.5 m <sup>3</sup> /h / 2bg	2 m <sup>3</sup> /h / 2bg	3 m <sup>3</sup> /h / 2bg

*Do you have some special requirements to fulfil?*

*We adapt ourselves to your process in order to provide you the tailor-made solution which suits best to your needs.*